

VIDAS® CAM

Campylobacter Enrichment Procedures



Carcass Rinse



Rinse with 100 ml of Bolton/Campylobacter Enrichment Broth with Laked Horse Blood and selective supplement for 1 minute

Incubate carcass rinse for 4 hours at 37°C in a microaerophilic atmosphere

Continue incubation for 20 hours at 42°C in a microaerophilic atmosphere

All Food Products



Add 25 g (25 ml) sample to 225 ml of Bolton/Campylobacter Enrichment Broth with Laked Horse Blood and selective supplement

Homogenized sample

Incubate 6 hours at 37±1°C in a microaerophilic atmosphere

Continue incubation for 42 hours at 42±1°C in a microaerophilic atmosphere

Store at 2-8°C for confirmation

1 - 2 ml



Heat at 95-100°C / 15±1 min

500 µl

VIDAS CAM



Results within 70 min

Store at 2-8°C for confirmation

VIDAS® CAM		
30111	VIDAS CAM	30 Tests
663183	Oxoid-Bolton Broth Selective Supplement	
R54072	Laked Horse Blood	
7526	Campylobacter Enrichment Broth (Base)	

For further details, consult the product package insert.

Test ⊖ < 0.10
Absence
of Campylobacter

Test ⊕ ≥ 0.10
CONFIRMATION

bioMérieux, Inc.
595 Anglum Road
Hazelwood, MO 63042
U.S.A.
Tel: (800) 634 7656
Fax: (800) 657 3053

bioMérieux Canada, Inc.
7815 Henri-Bourassa Boulevard West
St. Laurent (Québec)
Canada H4S 1P7
Tel: (1) 514 336 7321
Fax: (1) 514 807 0015

www.biomerieux-usa.com

