






GENERAL PROTOCOL



VIDAS UP Listeria (LPT) First Action Official Methods of Analysis #2013.10

Certified ISO 16140 by AFNOR (BIO 12/33-05/12) in human food products and production environment samples.

1. DETECTION USING VIDAS LPT FOR 125 G

 <p>Food Products</p> <p>125 g or 125 mL</p>	
<p>375 mL LPT Broth</p>  <p>24-30 hrs/ 30°C ± 1°C</p> <p>1 mL</p>	
 <p>1 ml into 10 ml LPT Broth</p> <p>22-26 hrs/ 30 ± 1 °C</p>	
<p>0.5 mL</p>  <p>VIDAS® Heat & Go / 5 ± 1* min</p>	
<p>Remove The Strip And Let It Cool For 10 min</p>	
 <p>VIDAS LPT Results Within 62 Minutes</p>	
<p>Test (-) <0.05 Absence of Listeria spp.</p>	<p>Test (+) ≥ 0.05 CONFIRMATION</p>

2. CONFIRMATION

LPT Broth		
Isolation on Selective Agar		
 <p>Regular blue/green colony. Presence of a typical halo for <i>Listeria monocytogenes</i>.</p> <p>ALOA® Agar 24 to 48 hrs/ 37°C ± 1°C</p>		
<p>VITEK 2GP, OMA 2012.02</p>	 <p>API® Listeria</p>	<p>Reference Method Confirmation</p>

LISTERIA RANGE

30126	VIDAS UP Listeria (LPT)	60 tests
410845	LPT Broth – Tubes	20 x 10 mL
410846	LPT Broth – Bottles	6 x 100 mL
410848	LPT Broth – Bottles	6 x 225 mL
410847	LPT Broth MiniBags	10 x 225 mL
410849	LPT Broth	4 x 3 L
AEB520080	ALOA® Agar 90 mm	20 plates
AEB520079	ALOA® Agar 90 mm	120 plates
AEB120082	ALOA® Agar 140 mm	10 plates
AEB15072	ALOA® Agar Dehydrated	500 g
M1091	Oxford(Mox) Agar	10 plates
10300	API® Listeria	10 strips
21342	VITEK 2GP	20 tests
93554	VIDAS® Heat & Go	110 volts

*For most of the samples, for further details, consult the product package insert.

**Consult us for availability in your country.