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# TRYPCASE SOY BROTH - Clean Rooms (TSB-SP)

*For microbiological control only*

Sterility testing by membrane filtration

## SUMMARY AND EXPLANATION

TRYPCASE SOY BROTH – Clean Rooms is recommended by the various pharmacopoeias for sterility testing. It is intended for the detection of aerobic bacteria. This medium complies with the European, US, and Japanese pharmacopoeias (1, 2, 3) for sterility testing by membrane filtration.

## PRINCIPLE

The broth contains a mixture of peptones which encourage the growth of most non-fastidious micro-organisms.

The bottles of TRYPCASE SOY BROTH – Clean Rooms are packaged in a plastic box which is double-wrapped and treated with ethylene oxide.

This packaging facilitates the introduction of the bottles into the protected area, limits the risk of contamination from the box, and avoids prior decontamination of the bottles.

An indicator enables visual confirmation that the ethylene oxide treatment has been performed.

Ethylene oxide treatment guarantees a 10<sup>6</sup> reduction in *Bacillus subtilis* spores on the outside of the bottles, the plastic box and the 2 bags.

This treatment does not modify the specifications of the medium.

## CONTENT OF THE KIT

	<b>Ready-to-use medium</b>
<b>REF 42 633</b>	6 x100 ml bottles in a plastic box Double-wrapped

## COMPOSITION

### Theoretical formula.

This medium can be adjusted and/or supplemented according to the performance criteria required:

Casein peptone (bovine).....	17 g
Soy peptone .....	3 g
Sodium chloride.....	5 g
Dipotassium phosphate .....	2.5 g
Dextrose.....	2.5 g
Purified water.....	1 l

pH 7.3

## REAGENTS AND MATERIAL REQUIRED BUT NOT PROVIDED

### Reagents:

- Fluid A (Ref. 42 044 or 42 616).
- Fluid D (Ref. 42 624).
- CLEAR THIO BROTH - Clean Rooms (Ref. 42 634)

### Material:

- Bacteriology incubator.

## WARNINGS AND PRECAUTIONS

- **For microbiological control only.**
- **For professional use only.**

- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest or inhale).
- All specimens, microbial cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling the bacterial group studied should be observed throughout this procedure. Refer to " CLSI® M29-A, *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline – Current Revision*". For additional information on handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories – CDC/NIH – Latest edition", or the current regulations in the country of use.
- Culture media should not be used as manufacturing material or components.
- Do not use reagents after the expiry date.
- Do not use bottles which show signs of contamination.
- Prior to use, decontaminate the outside of the bottles and the grey septum, if:
  - the lid of the plastic box does not include a round opening.
  - the bags show signs of damage.
  - a plastic cap covers the grey septum.
- According to the European Pharmacopoeia, the medium must be handled under Class A conditions (which should be situated within a class B environment) or in an isolator.
- The medium should be used according to the procedure indicated in this package insert. Any change or modification in the procedure may affect the results.

## STORAGE CONDITIONS

- **Store the bottles in their box at 2-25°C until the expiry date.**

## SPECIMENS

Follow the recommendations of the pharmacopoeias to perform specimen preparation.

## INSTRUCTIONS FOR USE

### Sterility testing by membrane filtration:

According to the method described in the pharmacopoeias, the following technique must be applied:

1. Filter a small quantity of Rinsing Solution (Fluid A or Fluid D) through each membrane.
2. Divide the test product into two equal parts and filter through separate membranes. The quantity of filter will depend on the type of product.
3. Rinse each membrane with 100 ml of Rinsing Solution. The number of rinses will depend on the product tested.
4. Pour 100 ml of TRYPCASE SOY BROTH – Clean Rooms onto one membrane and 100 ml of CLEAR - THIO BROTH - Clean rooms onto the other.
5. Incubate for a minimum of 14 days:
  - TRYPCASE SOY BROTH – Clean Rooms at **20-25°C**.
  - CLEAR - THIO BROTH - Clean Rooms at **30-35°C**.

### Direct inoculation:

Follow the methods described by the pharmacopoeias. It is possible to use the ready-to-use medium in tubes (Ref. 42 100) (see the corresponding package insert).

## READING AND INTERPRETATION

Refer to the indications of the pharmacopoeia. After incubation, observe the bacterial growth which is generally associated with turbidity of the broth.

## QUALITY CONTROL

TRYPCASE SOY BROTH – Clean Rooms is designed and developed to meet the strictest quality requirements. The results of the strains tested in the batch by batch quality control are given on the quality control certificate available on request.

## LIMITATIONS OF THE METHOD

- Growth depends on the requirements of each individual microorganism. It is therefore possible that certain strains which have specific requirements (substrate, temperature, incubation conditions, etc.) may not develop.
- Given the wide variety of specimens tested, it is the responsibility of the user to validate this broth for its specific intended use.
- As the color changes observed with ethylene oxide treatment indicators are variable, this result cannot be used as a criterion for batch release. A batch is released if the results of the biological indicators are conform.

## WASTE DISPOSAL

Dispose of used and unused reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products. It is the responsibility of each laboratory to handle waste and effluents produced according to their nature and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

## LITERATURE REFERENCES

1. European Pharmacopoeia 6.
2. USP 30.
3. Japanese Pharmacopoeia 15.


## INDEX OF SYMBOLS

Symbol	Meaning
	Catalogue number
	Manufacturer
	Temperature limitation
	Use by
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests

## WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

CLSI is a used, pending, and/or registered trademarks belonging to Clinical and Laboratory Standards Institute Inc.

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Printed in France

bioMérieux and the blue logo are used, pending, and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.

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