



FOOD
DECIDE WITH
SAFETY CONFIDENCE



VIDAS® *L. monocytogenes* Xpress (LMX) Ultra Performance Summary

Food safety professionals have a variety of test kits to choose from when looking to fill a need in their laboratory. One of the main criteria used in making their decision is an evaluation of the certifications the test kit has received. Validation studies provide a user with confidence in the performance of the test kit. This document summarizes the certifications granted to the bioMérieux VIDAS® *L. monocytogenes* Xpress (LMX) test method. The VIDAS LMX method was certified by the AOAC Research Institute *Performance Tested Method*SM (PTM) program after the completion of a single laboratory study and an independent laboratory study. These studies included an evaluation of the following parameters: inclusivity, exclusivity, robustness, lot-to-lot/stability and method comparison. The VIDAS LMX method was also adopted as an AOAC Official Method of Analysis (OMA), First Action status after the completion of a multi-laboratory method comparison study. Additionally, the VIDAS LMX method has been certified NF VALIDATION (AFNOR) as an alternative method for the detection of *Listeria monocytogenes* through the completion of an independent single laboratory study and a multi-laboratory study. All of the validation studies demonstrated that the VIDAS LMX method was statistically equivalent to the corresponding reference method for the matrices tested with 95% confidence.

AOAC Performance Tested Method: PTM#091103

In September 2011, the VIDAS LMX Assay was granted *Performance Tested Method*SM (PTM) status by the AOAC Research Institute for the detection of *Listeria monocytogenes* in a variety of foods. The VIDAS LMX method was validated according to harmonized PTM and *Official Methods of Analysis* guidelines (2002). Results of the validation study demonstrated the ability of the VIDAS LMX method to: 1) detect 60 different *L. monocytogenes* strains; 2) correctly show negative results for 31 non-*L. monocytogenes* strains; 3) support a 12-month shelf life and quality of the test kit; 4) perform appropriately after varied protocol parameters including enriched sample boiling time (4, 5, 6 min), sample temperature after boiling (10, 25, 50°C) and time reagents held at room temperature (after refrigeration) prior to performing VIDAS LMX method (0, 30, 60, 120 min); and 5) demonstrate no statistically significant differences when compared to a reference method for 7 claimed matrices using a 25g test portion size: deli ham, processed American cheese, vanilla ice cream, frozen spinach, cooked shrimp, smoked white fish, creamy peanut butter. Independent testing included method comparison studies for 3 matrices: deli ham, vanilla ice cream, cooked shrimp.

Through the AOAC Research Institute GovVal program, five matrices were evaluated for their performance compared to the USDA/FSIS MLG and the MFHPB-30 references methods: deli turkey, deli ham, liver pate, raw fermented sausage, turkey hot dogs. This study was validated according to a protocol written following the AOAC Microbiology Guidelines (2012) and approved by AOAC, Health Canada and Canadian Food Inspection Agency representatives.

In January 2013, the VIDAS LMX Assay matrix extension study was approved by the AOAC Research Institute when 4 matrices (deli ham, deli turkey, ground beef and queso fresco) were evaluated using a 125g test portion size. This study was validated according to a protocol written using the AOAC Microbiology Guidelines (2012).

AOAC Official Method: Certificate No: 2013.11

In June 2013, the VIDAS LMX Assay was granted First Action AOAC *Official Methods of Analysis*SM (OMA) status by the AOAC Research Institute for the detection of *L. monocytogenes* after the completion of the AOAC PTM evaluation and a collaboratively studied evaluation of raw queso fresco at both 25g and 125g. Fourteen different laboratories received twelve replicate test portions from each of the three contamination levels of matrix were analyzed by VIDAS LMX and the AOAC OMA 993.12 reference method. Statistical analysis was conducted according to the probability of detection model and showed no statistically significant difference in the number of positive test portions detected by the VIDAS LMX method and the OMA 993.12 method.

AFNOR – NF Validation: Certificate No: BIO 12/27 – 02/10

In April 2010, the VIDAS LMX Assay was certified NF VALIDATION as an alternative analysis method for the detection of *L. monocytogenes* in food products for humans and environmental samples. This validation was obtained by comparison with the reference method described in the international standard ISO 11290-1 according to the standard EN ISO 16140 (2003). In the single laboratory study, the results showed a relative sensitivity level of 90.5% for the VIDAS LMX method for all food and environmental surface testing. In the interlaboratory study, twelve different laboratories evaluated eight replicate test portions from each of the three contamination levels. Each test portion was analyzed by the VIDAS LMX and the ISO 11290-1 reference method. All presumptive positive test portions were confirmed following plating on selective agar and identification. The results showed a relative sensitivity level of 97.9% for the VIDAS LMX method.

Table of Contents

VIDAS® LMX	1
Table A. Validation Study Technical Requirements	3
AOAC PTM Validation Study	4
Table B. AOAC PTM Method Developer Method Comparison Study Details	6
Table C. AOAC PTM Method Comparison Study Results	7
Table D. Complete Inclusivity List	8
Table E. Complete Exclusivity List	9
Official Methods of Analysis (OMA) Validation Study	5
Table F. AOAC OMA Collaborative Study Data Summary	10
AFNOR NF Validation Study	11
Table G. AFNOR Independent Expert Laboratory Study Summary	12
Table H. AFNOR Independent Expert Laboratory Study Data Summary – Relative Sensitivity	13
Table I. AFNOR Independent Expert Laboratory Study Data Summary – Relative Detection Limit	13
Table J. AFNOR Inter-laboratory Study – Testing Results	14
Table K. AFNOR Inter-laboratory Study – Relative Sensitivity	14
Glossary of Terms	15

VIDAS® *L. monocytogenes* Xpress (LMX)

Catalog Number - REF 30 123



Contents of the VIDAS LMX Kit


- 60 LMX Strips
- 60 LMX SPR®s
- Standard LMX
- LMX Positive Control
- Negative Control
- 1 MLE Card (Master Lot Entry)
- 1 Package insert provided in the kit or downloadable from www.biomerieux.com/techlib

Principle of the Assay


VIDAS® *Listeria monocytogenes* Xpress (LMX) method is for use on the automated VIDAS instrument for the detection of *Listeria monocytogenes* antigens using the ELFA method (Enzyme Linked Fluorescent Assay). The Solid Phase Receptacle (SPR®) serves as the solid phase as well as the pipetting device. The interior of the SPR is coated with proteins specific for *Listeria monocytogenes* receptors. Reagents for the assay are ready-to-use and pre-dispensed in the sealed reagent strips. All assay steps are performed automatically by the instrument. The reaction medium is cycled in and out of the SPR several times. An aliquot of enrichment broth is dispensed into the reagent strip. The *Listeria monocytogenes* receptors present will bind to the interior of the SPR. Unbound components are eliminated during the washing steps. The proteins conjugated to the alkaline phosphatase are cycled in and out of the SPR and will bind to any *Listeria monocytogenes* receptors, which are bound to the SPR wall. A final wash step removes unbound conjugate. During the final detection step, the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR. The conjugate enzyme catalyzes the hydrolysis of the substrate into a fluorescent product (4-Methyl-umbelliferone) the fluorescence of which is measured at 450 nm. At the end of the assay, results are automatically analyzed by the instrument, which calculates a test value for each sample. This value is then compared to internal references (thresholds) and each result is interpreted as positive or negative.

Figure 1. Flow diagrams showing the VIDAS LMX standard protocols

VIDAS LMX Standard Protocol – 25g



VIDAS LMX Standard Protocol – 125g



Performance Tested MethodSM (PTM)

The AOAC Performance Tested MethodsSM (PTM) program was formed in 1992 and is a method certification program for proprietary methods. Methods certified as Performance TestedSM were found to perform according to the manufacturer's documented claims and are used throughout the global marketplace and within the regulatory arena. The PTM program offers certification as an endpoint for method evaluation or as an entry to method validation for programs requiring increased confidence and method reproducibility information. Validation study protocols are written according to AOAC Microbiology Guidelines (2002 or 2012) and include the following technical requirements: inclusivity/exclusivity, method developer method comparison, independent laboratory method comparison, robustness, product consistency, product stability and instrument variation (where applicable) studies. More information can be found at www.aoac.org.

Official Methods of AnalysisSM (OMA)

The Official Methods of AnalysisSM (OMA) program is offered by AOAC INTERNATIONAL and evaluates chemistry, microbiology, and molecular biology methods. It also evaluates traditional bench-top methods, instrumental methods, and proprietary, commercial, and/or alternative methods. Validation study protocols are written according to AOAC Microbiology Guidelines (2002 or 2012) and include the following technical requirements: inclusivity/exclusivity, method developer method comparison and collaborative method comparison studies. More information can be found at www.aoac.org.

AFNOR Validation Study

The NF VALIDATION mark has been widely recognized in France since the 1990s and is now well-established in Europe and internationally. It is a completely separate European certification system, operating alongside the technical validation systems of NordVal (inter-governmental validation system of 5 Nordic countries) and AOAC (North American technical validation system). Validation study protocols are written according to EN ISO 16140 and include the following technical requirements: inclusivity/exclusivity, single laboratory methods comparison and an inter-laboratory studies. More information can be found at <http://www.afnor-validation.com/afnor-validation-food-industry/food-industry.html>.

GovVal Study

The GovVal program was based on the PTM program operated by the AOAC Research Institute. Candidate test kit methods were evaluated using a validation protocol approved by AOAC, the Canadian Food Inspection Agency and Health Canada, using blind coded, randomized samples. Sample preparation and independent testing of many test kit methods, as well as the MFHPB-30 reference method, was done by an AOAC RI approved independent laboratory. The GovVal program was designed to evaluate previously AOAC-approved methods for the specific needs of regulatory agencies to enforce their standards for regulatory testing, in this case *Listeria monocytogenes* in ready-to-eat meats and *Listeria* species on stainless steel surfaces as compared to the Health Canada reference method.

Table A. Validation Study Technical Requirements

Study Type	AOAC PTM	AOAC OMA	AFNOR
Method Developer	●	●	
Independent	●		●
Collaborative		●	●

AOAC PTM Validation Study

PTM Certification#: 091103

PTM Certified: September 2011

PTM Matrix Extension Certified: January 2013

Guideline document – Original Study: AOAC INTERNATIONAL Methods Committee Guidelines for Validation of Qualitative and Quantitative Food Microbiological Official Methods of Analysis (2002)

Guideline document – Matrix Extension and GovVal: AOAC INTERNATIONAL Methods Committee Guidelines for Validation of Microbiological Methods for Food and Environmental Surfaces (2012)

Reference methods: USDA/FSIS Microbiological Laboratory Guidebook (8.07), US-FDA Bacteriological Analytical Manual Chapter 10 (2003), AOAC Official Method of Analysis 993.12 and Microbiology Food Health Protection Branch (MFHPB-30, Feb 2011)

Method Comparison (Tables B and C): The original method comparison study was performed at both the method developer (internal) and independent laboratories. Matrices (6 internal and 3 independent) were inoculated with *L. monocytogenes*. For each matrix, twenty replicates at one inoculation level (0.2-2 cfu/25g) and 5 uninoculated replicates were tested by both the VIDAS LMX and appropriate reference method. Primary enrichments for each method were confirmed using the traditional confirmation methods and ALOA chromogenic agar.

The data collected in the GovVal study was incorporated into the VIDAS LMX PTM study. There were 5 matrices evaluated using twenty replicates at two inoculation levels (0.2-2 cfu/test portion and 2-5 cfu/test portion) and 5 uninoculated replicates that were tested by the VIDAS LMX and the reference methods, USDA/FSIS MLG and MFHPB-30. Primary enrichments for each method were confirmed using the traditional confirmation methods and ALOA chromogenic agar followed by identification.

A matrix extension study was performed on four matrices at a larger test portion size, 125g. Each matrix was evaluated using twenty replicates at two inoculation levels (0.2-2 cfu/test portion and 2-5 cfu/test portion) and 5 uninoculated replicates that were tested by the VIDAS LMX and the appropriate reference method. Primary enrichments for each method were confirmed using the traditional confirmation methods and ALOA chromogenic agar followed by identification.

In the method comparison studies there were no significant differences between the VIDAS LMX method and the reference methods using unpaired Chi-square or the POD test at 5% level for the majority of the matrices evaluated. The only matrix that showed a significant difference, with the VIDAS LMX method resulting in a higher number of positive test portions, was turkey hot dogs (25g).

Inclusivity/Exclusivity (Tables D and E): The inclusivity studies demonstrated that the VIDAS LMX method could detect all 50 *L. monocytogenes* strains. For the exclusivity testing, all 30 non-*L. monocytogenes* strains tested negative by the VIDAS LMX method.

Lot-to-lot/Stability: Stability and lot-to-lot variation of the VIDAS LMX method was evaluated over 12 months at 2-8°C using *L. monocytogenes* and non-*L. monocytogenes* strains on three different test kit lots. There was no loss of stability over the 12 months or any difference between the three different test kits supporting the 12-month shelf life and the quality of the product.

Ruggedness: Minor variations to the protocol parameters, including enriched sample boiling time (4, 5, 6 min), sample temperature after boiling (10, 25, 50°C) and time reagents held at room temperature (after refrigeration) prior to performing VIDAS LMX method (0, 30, 60, 120 minutes). There were no differences seen in number of positive results after with each of the variations in the protocol.

Official Methods of Analysis (OMA) Validation Study

Certificate No: 2013.11

Certification date: June 2013

Guideline document: AOAC INTERNATIONAL Methods Committee Guidelines for Validation of Microbiological Methods for Food and Environmental Surfaces (2012)

Reference Method: AOAC 993.12 *Listeria monocytogenes* in Milk and Dairy Products

Inclusivity/Exclusivity & Method Developer Method Comparison

See the AOAC PTM study sections for evaluation results of these studies.

Collaborative Study: Queso fresco - 25g (Table F): Queso fresco (soft Mexican cheese, 25g) was analyzed by thirteen laboratories. The matrix was artificially contaminated with *L. monocytogenes* at two levels: a high level 3.81 CFU/25g (95% confidence interval of 3.06, 5.48) and a low level of 0.55 CFU/25g (95% confidence interval of 0.43, 0.70). A set of uncontaminated control test portions were also included at 0 CFU/25g. Twelve replicate test portions from each of the three contamination levels of matrix were analyzed. Two sets of test portions (72 total) were sent to each laboratory for analysis by VIDAS LMX and the AOAC 993.12 reference method due to different sample enrichments for each method. All test portions were confirmed by streaking the non-heated enrichment broth to ALOA and Oxford agar followed by identification.

Statistical analysis was conducted according to the probability of detection (POD). Results obtained by the different confirmation protocols showed no significant differences between the options evaluated.

Collaborative Study: Queso fresco - 125g (Table F): Queso fresco (soft Mexican cheese, 125g) was analyzed by twelve laboratories. The matrix was artificially contaminated with *Listeria monocytogenes* at two levels: a high level 5.41 CFU/125g (95% confidence interval of 3.53, 8.30) and a low level of 0.59 CFU/125g (95% confidence interval of 0.46, 0.74). A set of uncontaminated control test portions were also included at 0 CFU/25g. Twelve replicate test portions from each of the three contamination levels of matrix were analyzed. Two sets of test portions (72 total) were sent to each laboratory for analysis by VIDAS LMX and the AOAC 993.12 reference method due to different sample enrichments for each method. All test portions were confirmed by streaking the non-heated enrichment broth to ALOA and Oxford agar followed by identification.

Statistical analysis was conducted according to the probability of detection (POD). Results obtained by the different confirmation protocols showed no significant differences between the options evaluated.

Table B. AOAC PTM Method Developer Method Comparison Study Details

Matrix	Inoculating Organism	Test portion size(s)	Reference Method Comparison
Liver paté ^b	<i>L. monocytogenes</i> 1/2a, ATCC 49594 ^d	25g	USDA/FSIS MLG 8.07 & MFHPB-30
Raw fermented sausage ^b	<i>L. monocytogenes</i> 3b, CWD 1591 ^d	25g	USDA/FSIS MLG 8.07 & MFHPB-30
Turkey hot dogs ^b	<i>L. monocytogenes</i> 1/2b, ATCC BAA-751 plus 10x <i>L. innocua</i> 4ab ^d	25g	USDA/FSIS MLG 8.07 & MFHPB-30
Sliced deli turkey ^b	<i>L. monocytogenes</i> 3c, FSL-J1_049 ^d	25g	USDA/FSIS MLG 8.07 & MFHPB-30
Sliced deli turkey ^a	<i>L. monocytogenes</i> 4b, ATCC 51776 ^d	125g	USDA/FSIS MLG 8.07
Sliced deli ham ^c	<i>L. monocytogenes</i> 4b, ATCC 19115 ^d	25g	USDA/FSIS MLG 8.07 & MFHPB-30
Sliced deli ham ^a	<i>L. monocytogenes</i> ATCC 7644 ^d	125g	USDA/FSIS MLG 8.07
Ground beef ^a	<i>L. monocytogenes</i> 4b, ATCC 13932	125g	USDA/FSIS MLG 8.07
Queso fresco ^a	<i>L. monocytogenes</i> 4a, ATCC 19114	125g	AOAC OMA 993.12
Vanilla ice cream	<i>L. monocytogenes</i> 1	25g	AOAC OMA 993.12
Processed cheese	<i>L. monocytogenes</i> 1	25g	AOAC OMA 993.12
Frozen spinach	<i>L. monocytogenes</i> 4, Scott A	25g	FDA BAM Ch. 10 (2003)
Peanut butter	<i>L. monocytogenes</i> 4a	25g	FDA BAM Ch. 10 (2003)
Cooked shrimp	<i>L. monocytogenes</i> 3b	25g	FDA BAM Ch. 10 (2003)
Smoked white fish	<i>L. monocytogenes</i> 1/2b	25g	FDA BAM Ch. 10 (2003)

^aMatrix validated in PTM matrix extension study

^bMatrix validated as part of GovVal program

^cMatrix validated in both PTM and GovVal

^dHeat-stressed at 55°C for 10 minutes in a water bath

Table C. AOAC PTM Method Comparison Study Results

	INTERNAL LABORATORY DATA	Test portion size	Enrichment broth	Level	VIDAS LMX		Reference	χ²	Sensitivity, %	Specificity, %	False pos, %	False neg, %	dPOD CI	
					Presumptive	Confirmed							P vs C ^a	C vs R ^b
					Meat and Poultry	Deli ham ^a							25g	LMX
High	16	16	18	0.76			100	100	0	0	0 [-0.25, 0.25]	-0.10 [-0.33, 0.13]		
Deli ham	125g	LPT	Low	16		16	15	0.07	100	100	0	0	0 [-0.25, 0.25]	0.05 [-0.21, 0.29]
			High	15		15	17	0.61	100	100	0	0	0 [-0.26, 0.26]	0 [-0.26, 0.26]
Deli turkey ^a	25g	LMX	Low	15		15	15	0.00	100	100	0	0	0 [-0.26, 0.26]	0 [-0.26, 0.26]
			High	15		15	17	0.61	100	100	0	0	0 [-0.26, 0.26]	-0.1 [-0.34, 0.15]
Deli turkey	125g	LPT	Low	16		16	15	0.07	100	100	0	0	0 [-0.25, 0.25]	0.05 [-0.21, 0.29]
			High	16		16	15	0.07	100	100	0	0	0 [-0.25, 0.25]	0.05 [-0.21, 0.29]
Ground beef	125g	LPT	Low	7		7	7	0.00	100	100	0	0	0 [-0.24, 0.24]	0 [-0.24, 0.24]
			High	7		7	7	0.00	100	100	0	0	0 [-0.24, 0.24]	0 [-0.24, 0.24]
Liver paté ^a	25g	LMX	Low	11		11	10	0.10	100	100	0	0	0 [-0.28, 0.28]	0.05 [-0.24, 0.33]
			High	15		15	12	1.00	100	100	0	0	0 [-0.26, 0.26]	0.15 [-0.13, 0.40]
Raw fermented sausage ^a	25g	LMX	Low	9	9	10	0.10	100	100	0	0	0 [-0.28, 0.28]	-0.05 [-0.33, 0.24]	
			High	11	12	10	0.10	100	100	0	0	-0.05 [-0.33, 0.24]	0.05 [-0.24, 0.33]	
Turkey hot dogs ^a	25g	LMX	Low	10	10	8	0.39	100	100	0	0	0 [-0.28, 0.28]	0.10 [-0.19, 0.37]	
			High	20	20	13	8.27	100	100	0	0	0 [-0.16, 0.16]	0.35 [0.12, 0.57]	
Dairy	Processed American cheese	25g	LMX	Low	5	5	5	0.00	100	100	0	0	0 [-0.26, 0.26]	0 [-0.26, 0.26]
				High	5	5	5	0.00	100	100	0	0	0 [-0.26, 0.26]	0 [-0.26, 0.26]
	Vanilla ice cream	25g	LMX	Low	10	10	13	0.78	100	100	0	0	0 [-0.28, 0.28]	-0.15 [-0.41, 0.15]
Queso fresco	125g	LPT	Low	11	11	11	0.00	100	100	0	0	0 [-0.26, 0.26]	0 [-0.26, 0.26]	
Vegetable	Frozen spinach	25g	LMX	Low	4	4	5	0.21	100	100	0	0	0 [-0.25, 0.25]	-0.05 [-0.30, 0.21]
				High	4	4	5	0.21	100	100	0	0	0 [-0.25, 0.25]	-0.05 [-0.30, 0.21]
Seafood	Cooked shrimp	25g	LMX	Low	14	14	13	0.01	100	100	0	0	0 [-0.27, 0.27]	0.05 [-0.23, 0.32]
				High	14	14	13	0.01	100	100	0	0	0 [-0.27, 0.27]	0.05 [-0.23, 0.32]
Misc.	Smoked white fish	25g	LMX	Low	8	8	8	0.00	100	100	0	0	0 [-0.28, 0.28]	0 [-0.28, 0.28]
				High	8	8	8	0.00	100	100	0	0	0 [-0.28, 0.28]	0 [-0.28, 0.28]
Misc.	Creamy peanut butter	25g	LMX	Low	8	8	10	0.39	100	100	0	0	0 [-0.28, 0.28]	-0.10 [-0.37, 0.19]
				High	8	8	10	0.39	100	100	0	0	0 [-0.28, 0.28]	-0.10 [-0.37, 0.19]
	INDEPENDENT LABORATORY DATA	Test portion size	Enrichment broth	Level	VIDAS LMX		Reference	χ²	Sensitivity, %	Specificity, %	False pos, %	False neg, %	dPOD CI	
					Presumptive	Confirmed							P vs C ^a	C vs R ^b
						Deli ham							25g	LMX
High	16	16	14	0.52			100	100	0	0	0 [-0.25, 0.25]	0.10 [-0.17, 0.35]		
Low	8	8	12	1.56			100	100	0	0	0 [-0.28, 0.28]	-0.20 [-0.46, 0.10]		
	Vanilla ice cream	25g	LMX	Low	8	8	12	1.56	100	100	0	0	0 [-0.28, 0.28]	-0.20 [-0.46, 0.10]
				High	8	8	12	1.56	100	100	0	0	0 [-0.28, 0.28]	-0.20 [-0.46, 0.10]
	Cooked shrimp	25g	LMX	Low	7	7	6	0.11	100	100	0	0	0 [-0.28, 0.28]	0.05 [-0.23, 0.32]

^aTwo levels evaluated in the GoWal study

^aVIDAS LMX presumptive vs confirmed

^bVIDAS LMX vs reference method

$$\text{Mantel Haenszel Chi sq } (\chi^2) = \frac{(n-1)(ad-bc)^2}{(a+b)(a+c)(b+d)(c+d)}$$

N = total number of samples, a = candidate +, b = candidate -, c = reference +, d = reference -

Acceptability Criteria

χ² ≤ 3.84 indicates no significant difference (at the 0.05 level) between the two methods.

Sensitivity = VIDAS presumptive + (that confirmed +)/VIDAS confirmed +
Specificity = VIDAS presumptive - (that confirmed -)/VIDAS confirmed -
False positive = 100-sensitivity
False negative = 100-specificity

POD = x/N, where x is the number of positive test portions and N is the total number of test portions
dPOD = the difference between any two POD values

Acceptability Criteria

Confidence interval of a dPOD contains zero indicates no significant difference (at the 0.05 level) between the two methods

Table D. Complete Inclusivity List

	Strain	Origin		Strain	Origin
1	<i>Listeria monocytogenes</i> 1/2a	ATCC 35152	31	<i>Listeria monocytogenes</i> 1/2c	Ground beef
2	<i>Listeria monocytogenes</i> 1/2a	Smoked salmon	32	<i>Listeria monocytogenes</i> 1/2c	Beef in sauce
3	<i>Listeria monocytogenes</i> 1/2a	Pizza	33	<i>Listeria monocytogenes</i> 1/2	Smoked salmon
4	<i>Listeria monocytogenes</i> 1/2a	Raw milk cheese	34	<i>Listeria monocytogenes</i> 1/2	Chicken
5	<i>Listeria monocytogenes</i> 1/2a	Raw milk cheese	35	<i>Listeria monocytogenes</i> 3b	SLCC 2540
6	<i>Listeria monocytogenes</i> 1/2a	Rillettes (pate)	36	<i>Listeria monocytogenes</i> 3c	SLCC 2479
7	<i>Listeria monocytogenes</i> 1/2a	Raw milk cheese	37	<i>Listeria monocytogenes</i> 4a	ATCC 19114
8	<i>Listeria monocytogenes</i> 1/2a	Smoked salmon	38	<i>Listeria monocytogenes</i> 4b	Salad
9	<i>Listeria monocytogenes</i> 1/2a	Raw milk cheese	39	<i>Listeria monocytogenes</i> 4b	Raw milk cheese
10	<i>Listeria monocytogenes</i> 1/2a	Chicken breast	40	<i>Listeria monocytogenes</i> 4b	ATCC 19115
11	<i>Listeria monocytogenes</i> 1/2a	Ground beef	41	<i>Listeria monocytogenes</i> 4d	ATCC 19117
12	<i>Listeria monocytogenes</i> 1/2a	Sausage	42	<i>Listeria monocytogenes</i> 4e	ATCC 19118
13	<i>Listeria monocytogenes</i> 1/2a	Rabbit terrine	43	<i>Listeria monocytogenes</i> 4e	Reblochon (raw milk cheese)
14	<i>Listeria monocytogenes</i> 1/2a	Fried potatoes	44	<i>Listeria monocytogenes</i> 4e	Raw milk cheese
15	<i>Listeria monocytogenes</i> 1/2a	Fish	45	<i>Listeria monocytogenes</i> 7	SLCC 2482
16	<i>Listeria monocytogenes</i> 1/2a	Feed (soya)	46	<i>Listeria monocytogenes</i>	Sausage
17	<i>Listeria monocytogenes</i> 1/2a	Fried potatoes	47	<i>Listeria monocytogenes</i>	Salmon
18	<i>Listeria monocytogenes</i> 1/2b	Raw milk cheese	48	<i>Listeria monocytogenes</i>	Spinach
19	<i>Listeria monocytogenes</i> 1/2b	Pork tongue	49	<i>Listeria monocytogenes</i>	Cheese (Neufchatel)
20	<i>Listeria monocytogenes</i> 1/2b	Cream of chicken liver	50	<i>Listeria monocytogenes</i>	Mozzarella
21	<i>Listeria monocytogenes</i> 1/2b	Raw milk cheese	51	<i>Listeria monocytogenes</i>	Fish
22	<i>Listeria monocytogenes</i> 1/2b	SLCC 2755	52	<i>Listeria monocytogenes</i>	Vegetable
23	<i>Listeria monocytogenes</i> 1/2b	Pork ears	53	<i>Listeria monocytogenes</i>	Ham
24	<i>Listeria monocytogenes</i> 1/2c	Sausage from Montbéliard	54	<i>Listeria monocytogenes</i>	Ground beef
25	<i>Listeria monocytogenes</i> 1/2c	Ground beef	55	<i>Listeria monocytogenes</i>	Cheese (coulommier)
26	<i>Listeria monocytogenes</i> 1/2c	Beef	56	<i>Listeria monocytogenes</i>	Environment
27	<i>Listeria monocytogenes</i> 1/2c	Ground beef	57	<i>Listeria monocytogenes</i>	Environment
28	<i>Listeria monocytogenes</i> 1/2c	Pork belly	58	<i>Listeria monocytogenes</i>	Environmental sponges
29	<i>Listeria monocytogenes</i> 1/2c	Raw milk cheese	59	<i>Listeria monocytogenes</i>	French fries
30	<i>Listeria monocytogenes</i> 1/2c	Environmental sponge	60	<i>Listeria monocytogenes</i>	Beef

Table E. Complete Exclusivity List

	Strain	Origin
1	<i>Bacillus cereus</i>	Beets
2	<i>Bacillus mycoides</i>	Environmental sample
3	<i>Bacillus stearothermophilus</i>	Dairy
4	<i>Bacillus sphaericus</i>	Meat
5	<i>Streptococcus bovis</i>	Institut Pasteur Collection
6	<i>Enterococcus durans</i>	Meat
7	<i>Streptococcus equinus</i>	Institut Pasteur Collection
8	<i>Enterococcus faecalis</i>	ATCC 19433
9	<i>Enterococcus faecium</i>	Institut Pasteur Collection CIP 5433
10	<i>Streptococcus anginosus</i>	Institut Pasteur Collection
11	<i>Lactobacillus casei</i>	Dairy
12	<i>Jonesia denitrificans</i>	Institut Pasteur Collection
13	<i>Staphylococcus intermedius</i>	Institut Pasteur Collection
14	<i>Staphylococcus aureus</i>	Yogurt
15	<i>Staphylococcus epidermidis</i>	Yogurt
16	<i>Micrococcus</i>	Environmental sample
17	<i>Listeria seeligeri</i>	French fries
18	<i>Listeria grayi</i>	ATCC 25 401
19	<i>Listeria seeligeri</i>	Environmental sample
20	<i>Listeria welshimeri</i>	Salmon
21	<i>Listeria ivanovii</i>	Institut Pasteur Collection
22	<i>Listeria ivanovii</i>	Environmental sample
23	<i>Listeria ivanovii</i>	Institut Pasteur Collection
24	<i>Listeria innocua</i>	Raw milk cheese
25	<i>Listeria innocua 6b</i>	Ground beef
26	<i>Listeria ivanovii</i>	Institut Pasteur Collection
27	<i>Listeria seeligeri 1/2b</i>	Pork tongue
28	<i>Listeria welshimeri</i>	Sausage
29	<i>Listeria innocua</i>	Gorgonzola
30	<i>Listeria seeligeri</i>	Raw milk cheese
31	<i>Listeria ivanovii</i>	Ground beef

Table F. AOAC OMA Collaborative Study Data Summary

Matrix	Level (CFU/test portion)	VIDAS LMX		Reference (N=156)	χ ^{2a}	Sensitivity	Specificity	False pos	False neg	dLPOD ^b	
		Presumptive (N=156)	Confirmed (N=156)							P vs C ^a	C vs R ^b
Queso fresco (25g)	0.55 (0.43, 0.70)	77	75	69	0.23	97	100	3	0	0.01 [-0.10, 0.13]	0.04 [-0.08, 0.15]
	3.81 (3.06, 5.48)	156	156	153	1.50	100	100	0	0	0 [0.02, 0.02]	0.02 [-0.01, 0.06]
Matrix	Level (CFU/test portion)	VIDAS LMX		Confirmed (N=156)	χ ^{2a}	Sensitivity	Specificity	False pos	False neg	dLPOD ^b	
		Presumptive (N=156)	Confirmed (N=156)							P vs C ^a	C vs R ^b
Queso fresco (125g)	0.59 (0.46, 0.74)	70	70	69	0.01	100	100	0	0	0 [-0.12, 0.12]	0.10 [-0.10, 0.13]
	5.41 (3.53, 8.30)	144	144	144	0	100	100	0	0	0 [-0.03, 0.03]	0 [-0.03, 0.03]

^aVIDAS LMX presumptive vs confirmed

^bVIDAS LMX vs reference method

$$\text{Mantel Haenszel Chi sq } (\chi^2) = \frac{(n-1)(ad-bc)^2}{(a+b)(a+c)(b+d)(c+d)}$$

N = total number of samples, a = candidate +, b = candidate -, c = reference +, d = reference -

Acceptability Criteria

χ² ≤ 3.84 indicates no significant difference (at the 0.05 level) between the two methods.

Sensitivity = VIDAS presumptive + (that confirmed +)/VIDAS confirmed +

Specificity = VIDAS presumptive - (that confirmed -)/VIDAS confirmed -

False positive = 100-sensitivity

False negative = 100-specificity

POD = x/N, where x is the number of positive test portions and N is the total number of test portions

dPOD = the difference between any two POD values

Acceptability Criteria

Confidence interval of a dPOD contains zero indicates no significant difference (at the 0.05 level) between the two methods

AFNOR NF Validation Study

Certificate No: BIO 12/27 – 02/10

Certification date: February 2010

Guideline document: ISO 16140, Microbiology of food and animal feeding stuffs – Protocol for the validation of alternative methods (2003)

Reference Methods: ISO 11290-1, Microbiology of the food chain – Horizontal method for the detection and enumeration of *Listeria monocytogenes* and other *Listeria* spp.

Independent Expert Laboratory Study

Inclusivity/exclusivity (Tables D and E): The inclusivity studies demonstrated that the VIDAS LMX method could detect 50/50 *L. monocytogenes* strains tested after enrichment in LMX broth. For the exclusivity testing, 30/30 non-*L. monocytogenes* strains tested negative by the VIDAS LMX method.

Relative sensitivity (Tables G and H): The purpose of these tests was to evaluate the performance of the VIDAS LMX method with respect to the ISO 11290-1 reference method, on test portions naturally and artificially contaminated with *L. monocytogenes*, for the categories falling within the scope. The study evaluated 384 test portions that were inoculated with a wide variety of *L. monocytogenes* strains and analyzed at a level to achieve approximately 50% positive test portions with the following breakdown: meat – 62, dairy – 125, vegetables – 71, seafood – 66, environmental – 60. The results demonstrated that there are no statistically significant differences between the alternative method (90.5%) and the reference method (90%) when analyzing sensitivity values.

Relative level of detection (Table I): The objective of these tests was to determine the level of contamination needed to obtain about 50% of positive results and 50% negative results. Various “food matrix-strain” pairs were studied in parallel with the reference method and the VIDAS LMX method, for the studied categories. The results demonstrated that there are no statistically significant differences between the alternative method (0.2-1.8 cfu/25g) and the reference method (0.2-1.3 cfu/25g).

Inter-laboratory Study (Tables J and K): Pasteurized milk (25mL) was analyzed by twelve laboratories in this inter-laboratory study. The matrix was artificially contaminated with *L. monocytogenes* at two levels: a high level of 22.2 CFU/25g (95% confidence interval of 14, 33) and a low level of 3.26 CFU/25g (95% confidence interval of 0.7, 9.5). A set of uncontaminated control test portions were also included for each matrix at 0 CFU/25g. Eight replicate test portions from each of the three contamination levels of matrix were analyzed. Two sets of test portions (48 total) were sent to each laboratory for analysis by VIDAS LMX and the ISO 11290-1 reference method due to different test portion enrichments for each method. All test portions were confirmed following plating on selective agar and identification. The results demonstrated that there are no statistically significant differences between the alternative method (97.9%) and the reference method (100%) when analyzing sensitivity values.

Table G. AFNOR Independent Expert Laboratory Study Summary

	Matrix	Protocol Type	Enrichment broth	Test portion size	Number Positive	Number Negative	Total
Poultry	Raw	General	LMX	25g	9	11	20
	Raw seasoned	General	LMX	25g	10	10	20
	Delicatessen	General	LMX	25g	11	11	22
	TOTALS				30	32	62
Other meat	Raw milk & pasteurized milk chssese (cow)	General	LMX	25g	12	5	17
	Pasteurized milk cheeses (goat, ewe)	General	LMX	25g	13	13	26
	Raw milk cheese	General	LMX	25g	33	28	61
	Desserts, milk powders	General	LMX	25g	9	12	21
TOTALS				67	58	125	
Environment	Fish fillets & shellfish	General	LMX	25g	13	3	16
	Smoked fish	General	LMX	25g	10	24	34
	Ready-made meal	General	LMX	25g	8	8	16
TOTALS				31	35	66	
	Frozen	General	LMX	25g	9	12	21
	Raw/fresh	General	LMX	25g	13	13	26
	Seasoned	General	LMX	25g	10	14	24
	TOTALS				32	39	71
	Various waters	General	LMX	25g	9	10	19
	Surface samples	General	LMX	25g	11	9	20
	Residues & scraps	General	LMX	25g	10	11	21
TOTALS				30	30	60	
OVERALL TOTALS					190	192	384

Table H. AFNOR Independent Expert Laboratory Study Data Summary – Relative Sensitivity

	PA	NA	ND	PD	N	Relative Sensitivity, %	
						VIDAS LMX	Reference
All products	153	194	18	19	384	90.5	90
All foods (except raw milk cheese and environmental)	105	134	10	12	261	91.1	91.1
Totals:	258	328	28	31	645		
	PA	NA	ND	PD	N		
Meat products	24	32	3	3	62		
All Dairy	56	58	6	5	125		
Other dairy products	32	30	2	0	64		
Raw milk cheese	24	28	4	5	61		
Fish	27	35	2	2	66		
Vegetable	22	39	3	7	71		
Environmental	24	32	4	2	62		

PA – positive agreement (Candidate +/Reference +)
 NA – negative agreement (Candidate -/Reference -)
 ND – negative deviation (Candidate -/Reference +)
 PD – positive deviation (Candidate +/Reference -)
 N – total number of samples
 Relative Sensitivity, (PA + PD)/(PA + PD + ND)

Table I. AFNOR Independent Expert Laboratory Study Data Summary – Relative Detection Limit

Matrix	Test portion size	Strain	Relative detection level (CFU/25g or 25mL) with confidence interval ^a LOD ₅₀	
			Reference Method	VIDAS LMX Method
Rilette	25g	<i>L. monocytogenes</i> 1/2b	0.7 [0.4 - 1.3]	1.1 [0.6 - 1.8]
Raw milk	25g	<i>L. monocytogenes</i> 1/2b	0.5 [0.3 - 0.9]	0.4 [0.2 - 0.7]
Raw milk cheese	25g	<i>L. monocytogenes</i> 1/2b	0.7 [0.4 - 1.3]	0.6 [0.4 - 1.0]
Worn red cabbage	25g	<i>L. monocytogenes</i> 4b	0.5 [0.3 - 0.8]	0.6 [0.4 - 1.0]
Smoked salmon	25g	<i>L. monocytogenes</i> 1/2a	0.5 [0.3 - 0.9]	0.4 [0.2 - 0.7]
Process water	25g	<i>L. monocytogenes</i> 1/2c	0.4 [0.2 - 0.7]	0.3 [0.2 - 0.6]

^a LOD₅₀: estimated level of contamination enabling positive detection using the alternative method in 50% of cases.

Table J. AFNOR Inter-laboratory Study – Testing Results

Level	# positive/total	
	VIDAS LMX	Reference
Level 0 (0 CFU/25g)	0/96	0/96
Low level (3.26 CFU/25g)	96/96	96/96
High level (22.2 CFU/25g)	91/96	95/96

Table K. AFNOR Inter-laboratory Study – Relative Sensitivity

Level	PA	NA	ND	PD	N	Relative Sensitivity, %	
						VIDAS	Reference
Level 0 (0 CFU/25g)	-	96	-	-	96	97.9	100
Low level (3.26 CFU/25g)	96	-	-	2	96		
High level (22.2 CFU/25g)	92	-	4	1	96		

PA – positive agreement (Candidate +/Reference +)
 NA – negative agreement (Candidate -/Reference -)
 ND – negative deviation (Candidate -/Reference +)
 PD – positive deviation (Candidate +/Reference -)
 N – total number of samples
 Relative Sensitivity, (PA + PD)/(PA + PD + ND)

Glossary of Terms

Chi square – Test for significant difference; results less than 3.84 indicates no significant difference between methods.

Collaborative Study (AOAC) – A validation study performed by multiple laboratories to estimate critical VIDAS LMX method performance parameters. See also inter-laboratory study.

Confirmed result – The qualitative response from the confirmatory phase of the VIDAS LMX method.

False negative – A VIDAS negative test result that was confirmed to be culturally positive from the corresponding VIDAS enrichment (AOAC studies).

False positive – A VIDAS positive test result that was not confirmed culturally from the corresponding VIDAS enrichment (AOAC studies).

Fractional recovery – Validation criterion that is satisfied when an unknown sample yields both positive and negative responses within a set of replicate analyses. The proportion of positive responses should fall within 25% and 75% and should ideally approximate 50% of the total number of replicates in the set.

Inter-laboratory Study (AFNOR) – A validation study performed by multiple laboratories to estimate critical VIDAS LMX method performance parameters. See also collaborative study.

Limit of detection – The VIDAS system is able to detect a single target cell in a specific test portion following enrichment.

Negative agreement – Results for both the VIDAS LMX method and the reference method are negative.

Negative deviation – Result for the VIDAS LMX method is negative and the result for the reference method is positive.

Presumptive result – The qualitative response from the presumptive phase of the VIDAS LMX method that includes a confirmatory phase.

Probability of Detection (POD) – The proportion of positive analytical outcomes for a qualitative method for a given matrix at a given analyte level or concentration. POD is concentration dependent.

Positive agreement – Results for both the VIDAS LMX method and the reference method are positive.

Positive deviation – Result for the VIDAS LMX method is positive and the result for the reference method is negative.

Relative sensitivity – Ability of the alternative method to detect the analyte when it is detected by the reference method

Sensitivity – The number of VIDAS positive test results that were confirmed to be positive from the VIDAS enrichment divided by the total number of confirmed positive VIDAS enrichments.

Specificity – The number of VIDAS negative test results that were confirmed to be negative from the VIDAS enrichment divided by the total number of confirmed negative VIDAS enrichments.

Test portion – A specified quantity of the sample that is taken for analysis by the method.

For more information, please contact your local
bioMérieux Account Representative.
www.bioMerieux-usa.com/phage

