



VIDAS[®]

Immunoassay Reference Guide

OVER 30 YEARS OF PROVEN IMMUNOCHEMISTRY

The VIDAS family of compact immunoassay analyzers enhances automation, traceability, and productivity for your specialty testing.

- Random access anytime
- Ready-to-use reagents: 1 patient, 1 test, 1 result
- Robust systems: VIDAS 3 Mean Time Between Failure > 370 days

For more information, contact your local bioMérieux, Inc. representative.
Customer Service: (800) 682-2666 | www.biomerieux-usa.com/vidas-3

VIDAS® Assay Reference Sheet

ASSAY NAME	PRODUCT NUMBER	REPORTABLE RANGE	COMPATIBILITY*	CPT CODE (MODIFIER)	AVAILABLE ON	
					MINI VIDAS/VIDAS	VIDAS 3
CRITICAL AND EMERGENCY CARE						
VIDAS NEPHROCHECK® NEW!	421172-01	AKIRISK™ Score	N/A	N/A		✓
VIDAS B·R·A·H·M·S PCT™ (Procalcitonin)	30450-01	0.05-200 ng/ml	N/A	84145	✓	✓
VIDAS D-DIMER EXCLUSION™ II	30455-01	45-10,000 ng/ml	N/A	85379	✓	✓
INFECTIOUS DISEASES						
VIDAS MEASLES IgG	30219	Qualitative	1	86765	✓	✓
VIDAS MUMPS IgG	30218	Qualitative	1	86735	✓	✓
VIDAS RUBELLA IgG	30226	0-250 IU/ml	2	86762	✓	✓
VIDAS VARICELLA IgG	30217	Qualitative	1	86787	✓	✓
VIDAS LYME IgG II	417401	Qualitative	3	86618	✓	✓
VIDAS LYME IgM II	416436	Qualitative	3	86618	✓	✓
VIDAS <i>C. difficile</i> TOXIN A&B (CDAB)	30118-01	Qualitative	N/A	87324	✓	✓
VIDAS <i>C. difficile</i> GDH	30125-01	Qualitative	N/A	87449	✓	✓
VIDAS <i>H. pylori</i>	30192-01	Qualitative	2	86677	✓	✓
VIDAS TOXO IgG II	30210-01	0-300 IU/ml	2	86777	✓	✓
VIDAS TOXO IgM	30202-01	Qualitative	2	86778	✓	✓
VIDAS TOXO Avidity IgG	30222-01	Avidity Index	2	86777	✓	
VIDAS CMV IgG	30204-01	0-400 AU/ml	2	86644	✓	✓
VIDAS CMV IgM	30205-01	Qualitative	N/A	86645	✓	✓
HORMONES (REPRODUCTIVE AND THYROID)						
VIDAS HCG	30405-01	2-1500 mIU/ml	N/A	84702	✓	✓
VIDAS LH	30406-01	0.1-100 mIU/ml	4	83002	✓	
VIDAS FSH	30407-01	0.1-110 mIU/ml	4	83001	✓	
VIDAS PROGESTERONE	30409-01	0.25-80 ng/ml	N/A	84144	✓	
VIDAS ESTRADIOL II	30431-01	9-3000 pg/ml	N/A	82670	✓	
VIDAS TSH	30400-01	0.05-60 µIU/ml	N/A	84443	✓	

*Matching numbers indicate VIDAS assays that can be run using the same protocol on the VIDAS instrument.

To place an order, visit

www.biomerieuxDIRECT.com