



Variety of Foods:

The **VIDAS® *Listeria monocytogenes* Express (LMX) Assay #091103** has been validated and certified as a *Performance Tested Method[®]* by the AOAC Research Institute as an effective method for the detection of ***Listeria monocytogenes*** in a variety of foods. Test portion size validated was 25g.

Vegetables: spinach, peanut butter (BAM)

Seafood: shrimp, smoked fish (BAM)

Meat: liver paté, raw fermented sausage (USDA & MFHPB30)

Poultry: turkey hot dogs, deli ham, deli turkey (USDA & MFHPB30)

Dairy: processed cheese, ice cream (AOAC)