


# VIDAS® Easy SLM

## Detection and confirmation of Salmonella in food products and animal feed


Performance Tested™ by AOAC Research Institute (Certificate # 020901)

General protocol, validated ISO 16140 by AFNOR, N°BIO 12/16-09/05

### 1 detection using VIDAS Easy SLM



### 2 confirmation with reference method



#### Salmonella range

42121	Salmonella Xpress 2 broth	20 tubes 10 ml
43621	ChromID Salmonella	10 plates
M1031	XLD - agar	10 plates
M1016	Hektoen - agar	20 plates
20100	API 20 E	25 strips
32700	ID 32 E	25 strips
30702	VIDAS Salmonella (SLM)	60 tests
93554	VIDAS Heat and Go*	110 Volts

For further details, consult the product package insert.

\* Please contact your local representative.

bioMérieux, Inc.  
595 Anglum Road  
Hazelwood, MO 63042  
U.S.A.  
Tel: (800) 634 7656  
Fax: (800) 657 3053

bioMérieux Canada, Inc.  
7815 Henri-Bourassa Boulevard West  
St. Laurent (Québec)  
Canada H4S 1P7  
Tel: (1) 514 336 7321  
Fax: (1) 514 807 0015