


VIDAS® SLM

Detection and Confirmation of *Salmonella* (TT & RV)


Protocol validated by AOAC International (Official Method of Analysis, Method No. 2004.03)



1. Detection using VIDAS® SLM



2. Confirmation



VIDAS® SLM		
T8034	Tetrathionate Broth Base	20 x 10 ml
42073	Rappaport Vassiliadis Broth	20 x 10 ml
T8080	M Broth	20 x 10 ml
30702	VIDAS SLM	60 tests
M1016	Hektoen Enteric Agar	10 plates
M1031	XLD	10 plates
7113	Bismuth Sulfite Agar	500g
43621	ChromID™ Salmonella	20 plates
21341	VITEK 2 GN	20 cards
20100	api 20E	25 strips
32700	rapid ID 32 E	25 strips

For further details, consult the product package insert.

bioMérieux, Inc.
595 Anglum Road
Hazelwood, MO 63042
U.S.A.
Tel: (800) 634 7656
Fax: (800) 657 3053

bioMérieux Canada, Inc.
7815 Henri-Bourassa Boulevard West
St. Laurent (Québec)
Canada H4S 1P7
Tel: (1) 514 336 7321
Fax: (1) 514 807 0015

www.biomerieux-usa.com

