


**Safety & Quality  
Solutions**



**BacT/ALERT 3D**

**VITEK 2**  
—technology

**OBSERVA** 

**Integrated Data Management Software  
Connecting Your World**

  
**BIOMÉRIEUX  
INDUSTRY**



## Integrated Data Management

*A data management system designed specifically for the most demanding environment... Your laboratory!*

Every lab is unique – the sample types, workload, staffing, space, information system needs. The right microbiology solution for someone else’s laboratory might not be right for yours. Our design approach allows connectivity with either the BacT/ALERT® 3D microbial detection system or the VITEK® 2 Compact automated identification system – or both. OBSERVA® lets you configure your system the way you want, depending on your laboratory needs.



### Features

- Intuitive user interface
- Automatic data back up, archive to DVD
- Bi-directional interface
- Password protection
- Detailed audit trails
- Multiple access levels

### Specifications

- Windows® XPE, poet database
- HP d530 computer
- HP LIS02 LCD monitor
- HP 1320 laser printer
- UPS
- External modem

### Reports


#### BacT/ALERT 3D

- Status reports
- Load report
- Unload report
- Days to detection detail report
- Incidence detail report
- Statistical reports


#### VITEK 2 Compact

- Product report
- Accession report
- Organism trending report
- Bionumber distribution report
- System log event report

Password level screen shot



Reports screen shot



**bioMérieux, Inc.**  
 595 Anglum Road  
 Hazelwood, MO 63042-2320  
**U.S.A.**  
**Tel: (800) 634 7656**  
**Fax: (800) 657 3053**

**bioMérieux Canada, Inc.**  
 7815 Henri-Bourassa Boulevard West  
 St. Laurent (Québec)  
 Canada H4S 1P7  
**Tel: (1) 514 336 7321**  
**Fax: (1) 514 807 0015**

[www.biomerieux-usa.com](http://www.biomerieux-usa.com)

